


# Important advice

Prior to using the magnetic optiCLIC sticks of the GAMMAT® optimo2 or of the GAMMAT® free, it is absolutely necessary to carefully clean the thin end of the stick with a wire brush, grinding paper, or rubber polishers in order to remove the (invisible) oxide layer developing on the stainless steel.



⚠ Please make sure that the plastic coating is not harmed.



**Gramm Technik GmbH**  
**Dental Division**

Einsteinstrasse 4  
D-71254 Ditzingen-Heimerdingen  
Germany  
Phone +49 (0)7152 500960  
Fax +49 (0)7152 500980  
E-Mail [gramm.dental@t-online.de](mailto:gramm.dental@t-online.de)  
[www.electroforming.de](http://www.electroforming.de)

Revised: March 2008